

*Innovative Technology
for Non-Invasive
Ventilation*

BiTrac MaxShield™

Pulmonary® BiTrac MaxShield™ is ideal for use in Non-Invasive Ventilation (NIV), CPAP and pressure support Bilevel and Ventilator Applications.

BiTrac MaxShield's contoured cushion provides a quick and effective seal on the perimeter of the face.


The **BiTrac MaxShield** provides a perfect solution for patients who can not tolerate full face or nasal masks and is ideal for patients requiring immediate ventilation with a quick and simple set up by reducing the time of the fitting process.

Convenient *single use* eliminates cross contamination and issues relating to cleaning and disinfection.

- *BiTrac MaxShield* for CPAP, Bilevel and Ventilator Applications
- Silicone cushion for full facial perimeter low pressure seal
- Prevents nasal skin breakdown
- Unobstructed view for claustrophobic patients
- Two sizes fit most adults, simplifies mask fitting
- MRI Compatible
- **DISPOSABLE** - Single Patient Use
-  **Latex Free**
- *Made in USA*



*Complete line
of masks for
immediate
non-invasive
ventilation*



BiTrac MaxShield™



BiTrac SE MaxShield

BiTrac MaxShield

BiTrac MaxShield (Bilevel and CPAP Application)

with Anti-Asphyxia Elbow 22mm Female and MaxShield Head Strap
313-9052 Adult
313-9054 Small Adult

10 per case
10 per case

BiTrac SE MaxShield (Ventilator Application)

with Standard Elbow 22mm Female and MaxShield Head Strap
313-9053 Adult
313-9055 Small Adult

10 per case
10 per case

Circuits

BiTrac™NIV Circuit (Bilevel/CPAP Applications)

313-9074 72" Smooth ID Tubing with swivel connector with exhalation port, cap plug; includes optional 24" tubing with "T"

10 per case

BiTrac NIV Circuit (Bilevel/CPAP Applications)

313-9075 72" Smooth ID Tubing with attached filter, swivel connector with exhalation port, cap plug; includes optional 24" tubing with "T"

10 per case

ValueTrac™NIV Circuit (Bilevel/CPAP Applications)

313-9078 72" Corrugated Tubing and swivel connector with exhalation port

10 per case



BiTrac MaxShield

BiTrac MaxShield allows unobstructed view minimizing claustrophobia

Convenient, pre-attached Head Strap

Two exhaust ports provide continuous leak (BiTrac SE MaxShield does not have exhaust ports)

Manometer port for connecting proximal pressure line or for monitoring CO₂

**Bilevel
Compatible™**

Anti-Asphyxia Elbow for patient safety during Bilevel Applications (BiTrac SE MaxShield has a Standard Elbow for ventilator application)

Head Strap with quick and easy attachment clips



... bringing change to life®

Pulmonary
2055 Executive Drive
Indianapolis, IN 46241 USA

T: 317.246.5505
F: 317.246.5501

www.Pulmonary.com



0344

Pulmonary and bringing change to life are registered trademarks of Pulmonary, Inc. BiTrac, BiTrac MaxShield, Bilevel Compatible and ValueTrac are trademarks of Pulmonary, Inc.